ANICE COVID-19 Antigen Saliva Test Kit (Collaidal Gold) Ujian Kendiri

Arahan Penggunaan

ISP-CA1

[Nama produk]

ANICE COVID-19 Antigen Saliva Test Kit (Colloidal Gold)Ujian Kendiri

VIDEO TUTORIAL

[Tujuan Penggunaan]

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Ketepatan: 98.64% (97.57%–99.32%)

[Prestasi Analisis]

Had Pengesanan (Sensitiviti Analitik)

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ANICE COVID-19 Antigen Saliva Test Kit (Colloidal Gold)(Self-Test Only)

Instruction for Use

VIDEO TUTORIAL

REF ISP-CA1

【Product name】

ANICE COVID-19 Antigen Saliva Test Kit (Colloidal Gold)(Self-Test Only)

【Intended use】

COVID-19 Antigen Rapid Saliva Test Kit (Colloidal Gold)(Self-Test Only) is used for in vitro qualitative detection of the SARS-CoV-2 antigen (Nucleocapsid protein) in human saliva specimen. The novel corona virus belongs to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel corona virus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

【Test principle】

This kit uses immunochromatography for detection. The specimen will move forward along the test card under capillary action. If the specimen contains SARS-CoV-2 antigen, the antigen will bind to the colloidal gold-labeled new corona virus monoclonal antibody. The immune complex will be captured by corona virus monoclonal antibodies which are membrane fixed, form the fuchsia line in detection line, display will be SARS-CoV-2 antigen positive; if the line does not show color, and it means the negative result. The test card also contains a quality control line C, which shall appear fuchsia regardless of whether there is a detection line. Specifications and Main Components

【Specifications and Main Components】

- 1 COVID-19 Antigen Saliva Test Kit
- 1 Extraction Tube
- 1 Saliva Funnel
- 1 Disposable dropper
- 1 Instructions For Use
- 1 Biohazard disposal bag and 1 disinfectant wipes

Following materials required but not provided in the kit:

- Timer (watch)
- Any necessary personal protective equipment (gloves, glasses etc.)


【Storage and Shelf life】

1. Store at 2°C-30°C, and the shelf life is 18 months.
2. After the aluminum foil bag is unsealed, the test kit should be used as soon as possible within one hour.

【Specimen Collection】

Note: it is recommended to collect saliva samples after getting up in the morning. Do not eat, drink, smoke or chew gum 30 minutes before saliva collection.

1. Collect saliva samples with a saliva collector.

2. Specimen storage

It is recommended to test within 2 hours after collection. It can be stored for up to 2 months if frozen below -20°C.

【Test Methods】

The test method was colloidal gold. Please read the manual and the instrument operation manual carefully before use.

1. Open the package and take out the covid-19 antigen saliva test kit.
2. Place the extraction tube include collected saliva into the Tube Holder of the carton.
3. Open the lid and draw a tube of liquid with a disposable dropper. Drop 2 drops into the sample well of the test kit and start the timer.
4. Read the results within 20 minutes. Strong positive results can be reported within 20 minutes, however, negative results must be reported after 20 minutes, and the results after 30 minutes are no longer valid.

【Waste Disposal after Test Procedures】

1. Place all used devices into biohazard bag and seal biohazard bag
2. Dispose all used devices into normal household waste container in compliance with the applicable local regulations.
3. Wash or sanitize your hands again.

【Result Interpretation】

Negative result: if there is only a quality control line C, the detection line is colorless, indicating that SARS-CoV-2 antigen has not been detected and the result is negative. Negative result indicates that the content of the SARS-CoV-2 antigen in the sample is below the limit of detection or no antigen. Negative results should be treated as presumptive, and do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's recent exposures, history, and the presence of clinical signs and symptoms consistent with COVID-19, and confirmed with a molecular assay, if necessary, for patient management.


Positive result: if both the quality control line C and the detection line appear, SARS-CoV-2 antigen has been detected and the result is positive for antigen.

The positive results indicate the existence of SARS-CoV-2 antigen. It should be further diagnosed by combining the patient's history and other diagnostic information. The Positive results do not rule out bacterial infection or co-infection with other viruses. Pathogens detected are not necessarily the main cause of disease symptoms.

Caution: There is currently a suspicion of a COVID-19 infection. You are therefore encouraged to: To contact a doctor/general practitioner or the local health department immediately. You need to self-report your result to My Sejahtera https://tmysejahtera.malaysia.gov.my/help_en/selfreport/

- **Comply with local guidelines for self-isolation.**
- **Have a PCR confirmatory test performed.**

Invalid result: if the quality control line C is not observed, it will be invalid regardless of whether there is detection line (as shown in the figure below), and the test shall be conducted again. The invalid result indicates that the procedure is not correct or that the test kit is out of date or invalid. In this case, the package insert should be read carefully and repeat the test with a new test device. If the problem persists, discontinue using the test kit of this Lot number immediately and contact your local distributor.

【Limitations】

1. This reagent is only used for in vitro diagnosis.
2. This reagent is only used to detect human saliva specimen. The results of other specimens may be wrong.
3. This reagent is only used for qualitative detection and can't detect the level of novel corona virus antigen in the specimen.
4. This reagent is only a clinical auxiliary diagnostic tool. If the result is positive, it is recommended to use other methods for further examination in time and the doctor's diagnosis shall prevail.
5. If the test result is negative and clinical symptoms persist. It is recommended to repeat sampling or use other testing methods for testing. A negative result cannot preclude the possibility of exposure to or infection with SARS-CoV-2 virus at any time.
6. The test results of the test kits are for clinicians' reference only, and should not be used as the only basis for clinical diagnosis. The clinical management of patients should be comprehensively considered in combination with their symptoms/signs, medical history, other laboratory tests and treatment responses, etc.
7. Due to the limitation of the detection reagent methodology, the limit of detection of this reagent is generally lower than that of nucleic acid reagents. Therefore, the test personnel should pay more attention to the negative results and need to combine other test results to make a comprehensive judgment. It is recommended to use nucleic acid testing or virus isolation and culture identification methods to review negative results which have doubts.
8. Positive test results do not exclude co-infection with other pathogens.
9. False negative results may occur when the SARS-CoV-2 antigen level in the sample is lower than the detection limit of the kit or the specimen collection and transportation are not appropriate. Therefore, even if the test results are negative, the possibility of SARS-CoV-2 infection can't be ruled out.
10. Positive and negative predictive values are highly dependent on prevalence rates. Positive test results are more likely to represent false positive results during periods of little/no SARS-CoV-2 activity when disease prevalence is low. False negative test results are more likely when prevalence of disease caused by SARS-CoV-2 is high.
11. Analysis of the possibility of false negative results:
 - (1) Unreasonable specimen collection, transportation and processing, low virus titer in the sample, no fresh sample or freezing and thawing cycling of the sample may lead to false negative results.
 - (2) The mutation of viral gene may lead to changes in antigenic determinants, which lead to negative results.
 - (3) The research on the SARS-CoV-2 has not been completely thorough; the virus may mutate and cause the differences for best sampling time (virus titer peak) and sampling location. Therefore, for the same patient, we can collect samples from multiple locations or follow up for multiple times reduce the possibility of false negative results.
12. Monoclonal antibodies may fail to detect, or detect with less sensitivity, SARS-CoV-2 viruses that have undergone minor amino acid changes in the target epitope region.

【Warnings and Precautions】

1. Read the instructions carefully before using the kit, and strictly control the reaction time. If you do not follow the instructions, you will get inaccurate results.
2. Guard against moisture, do not open the aluminum platinum bag before it is ready for testing. Do not use the aluminum foil bag when it is damaged or the test card is damp.
3. Please use it within the validity period.
4. Balance all reagents and specimens to room temperature before use.
5. Do not replace the components in this kit with components in other kits.
6. Do not dilute the specimen for testing, otherwise you may get inaccurate results.
7. The kit shall be stored in strict accordance with the conditions specified in this manual. Please do not store the kit under freezing conditions.
8. The test methods and results must be interpreted in strict accordance with this specification.
9. Negative results will occur with this kit if the novel corona virus antigen in the specimen falls below the minimum detection limit for this kit.

【Clinical Performance】

A total of 807 Saliva Sample (395 COVID-19 PCR positive patient samples and 412 COVID-19 PCR negative patient samples) from individual symptomatic patients were collected according to proper procedure.

The performance of ANICE COVID-19 Antigen Saliva Test Kit compared to the PCR comparator method is shown in the table below.

Method	PCR Test		Total	
	Result	Positive		Negative
ANICE COVID-19 Antigen Saliva Test Kit (Colloidal Gold) (For Self-Test Only)	Positive	393	2	395
	Negative	9	403	412
	Total	402	405	807

Sensitivity: 97.76% (95.79%~98.97%)

Specificity: 99.51% (98.23%~99.94%)

Accuracy: 98.64% (97.57%~99.32%)

【Analytical Performance】

Limit of Detection (Analytical Sensitivity)

The limit of detection (LoD) for the COVID-19 Antigen Saliva Test kit was established using heat-inactivated SARS-CoV-2 virus and SARS-CoV-2 recombinant nucleocapsid protein. The confirmed LoD for direct swab was 1.7x10² TCID₅₀/mL with virus and 1ng/mL with recombinant nucleocapsid protein.

Cross Reactivity (Analytical Specificity)

Cross reactivity and potential interference of COVID-19 Antigen Saliva Test kit was evaluated. No cross-reactivity or interference was observed in the table below.

Pathogen	Concentration
HKU1	1.0×10 ⁵ pfu/mL
MERS	1.0×10 ⁵ pfu/mL
OC43	1.0×10 ⁵ pfu/mL
NL63	1.0×10 ⁵ pfu/mL
229E	1.0×10 ⁵ pfu/mL
Influenza A virusH1N1	1.0×10 ⁵ pfu/mL
Influenza A virusH3N2	1.0×10 ⁵ pfu/mL
Influenza A virusH5N1	1.0×10 ⁵ pfu/mL
Influenza A virusH7N9	1.0×10 ⁵ pfu/mL
Influenza B virus	1.0×10 ⁵ pfu/mL
Rhinovirus A	1.0×10 ⁵ pfu/mL
Rhinovirus B	1.0×10 ⁵ pfu/mL
Rhinovirus C	1.0×10 ⁶ pfu/mL
adenovirus 1	1.0×10 ⁵ pfu/mL
adenovirus 2	1.0×10 ⁵ pfu/mL
adenovirus 3	1.0×10 ⁵ pfu/mL
adenovirus 4	1.0×10 ⁵ pfu/mL
adenovirus 5	1.0×10 ⁵ pfu/mL
adenovirus 7	1.0×10 ⁵ pfu/mL
adenovirus 55	1.0×10 ⁵ pfu/mL
Enterovirus A	1.0×10 ⁵ pfu/mL
Enterovirus B	1.0×10 ⁵ pfu/mL
Enterovirus C	1.0×10 ⁵ pfu/mL
Enterovirus D	1.0×10 ⁵ pfu/mL
EB virus	1.0×10 ⁵ pfu/mL
measles virus	1.0×10 ⁵ pfu/mL
Human cytomegalovirus	1.0×10 ⁵ pfu/mL
Rotavirus	1.0×10 ⁵ pfu/mL
Norovirus	1.0×10 ⁵ pfu/mL
Mumps virus	1.0×10 ⁵ pfu/mL
Varicella zoster virus	1.0×10 ⁵ pfu/mL
Respiratory syncytial virus	1.0×10 ⁵ pfu/mL
Mycoplasma pneumoniae	1.0×10 ⁶ cfu/mL
Streptococcus pneumoniae	1.0×10 ⁶ cfu/mL
Streptococcus pyogenes	1.0×10 ⁶ cfu/mL
Staphylococcus epidermidis	1.0×10 ⁶ cfu/mL
Staphylococcus aureus	1.0×10 ⁶ cfu/mL
Pseudomonas aeruginosa	1.0×10 ⁶ cfu/mL
Escherichia coli	1.0×10 ⁶ cfu/mL
Bordetella pertussis	1.0×10 ⁶ cfu/mL
Legionella pneumophila	1.0×10 ⁶ cfu/mL
Haemophilus influenzae	1.0×10 ⁶ cfu/mL
Candida albicans	1.0×10 ⁶ cfu/mL

Interference Study

The following substances were designed to be added to the negative and weak positive samples to test the influence on the test results (see the table below). And no interference was observed in the following substances.

Material	Concentration
Whole Blood	4%(v/v)
Mucin	0.5%(v/v)
Sucrets (Dyclonin)	1.5mg/mL
Sucrets (Menthol)	1.5mg/mL
Chloraseptic (Menthol)	1.5mg/mL
Chloraseptic (Benzocaine)	1.5mg/mL
CVS Nasal Drops (Phenylephrine)	15%(v/v)
Afrin (Oxymetazoline)	15%(v/v)
CVS Nasal Spray (Cromolyn)	15%(v/v)
Nasal Gel (Oxymetazoline)	10%(v/v)
Zicam	5%(v/v)
Ricola (Menthol)	1.5mg/mL

Homeopathic (Alkalol)	1:10
Fisherman's Friend	1.5mg/mL
Sore Throat Phenol Spray	15%(v/v)
Tobramycin	4µg/mL
Mupirocin	10mg/mL
Fluticasone Propionate	5%(v/v)
Tamiflu (Osetamivir Phosphate)	5mg/mL

However, high concentration of hemoglobin has an effect on the background elimination, suggesting that hemolytic samples have an impact on the background fading of products, so it can affect the observation of weak positive samples.

Hook Effect

No significant hook or front band effect was shown when tested at the concentration of 3.4x10⁵TCID₅₀/mL of the heat-inactivated SARS-CoV-2 Virus samples and at the concentration of 1µg/mL of the recombinant nucleocapsid protein.

【Symbol】

Symbol	Description	Symbol	Description
	CE mark		In vitro diagnostic medical device
	Lot number		Date of manufacture
	Expiry date		Do not reuse
	Caution		Consult instructions for use
	European authorized representative		Manufacturer
	Stored at 2-30°C		Keep away from sunlight
	Keep dry		Don't use the product when the package is damaged
	Sterilized using ethylene oxide		

ISP WELLNESS SDN BHD (1433750-M)

No 47, Jalan Chung Ah Ming, 31650 Ipoh, Perak, Malaysia.

Tel: +(60)5-321 7673