## ANICE COVID-19 Antigen Saliva Lollipop (Ujian sendiri)

### Arahan Penggunaan

ISP-CA2

#### 【Tujuan Penguanaan】

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Nama virus	Tahap Ujian
Endemic human coronavirusNL63	≤9.87x10sPFU/mL
MERS Coronavirus	≤793PFU/mL
Influenza A virus (H3N2)	<8.82x10s PFU/mL
Influenza B virus (Malaysia )	≤2.92x10sPFU/mL
Rhinovirus	≤4.17x10sPFU/mL
Endemic human coronavirus OC43	≤1.0x10sPFU/mL
Endemic human coronavirus 229E	≤1.0x10sPFU/mL
Adenovirus	≤1.0x10sPFU/mL
Metapneumovirus	≤1.0x10sPFU/mL
Parainfluenza virus type 1	≤1.0x10sPFU/mL
Parainfluenza virus type 2	≤1.0x10sPFU/mL
Parainfluenza virus type 3	≤1.0x10sPFU/mL
Parainfluenza virus type 4	≤1.0x10sPFU/mL
Influenza A virus (H1N1)	≤1.0x10sPFU/mL
Enterovirus	≤1.0x10sPFU/mL
Respiratory syncytial virus respectively	≤1.0x10sPFU/mL
Haemophilus influenzae	≤1.0x10sCFU/mL
Streptococcus pneumoniae	≤1.0x10sCFU/mL
Streptococcus pyogenes	≤1.0x10sCFU/mL
Candida albicans	≤1.0x10sCFU/mL
Bordetella parapertussis	≤1.0x10sCFU/mL
Mycoplasma pneumoniae	≤1.0x10sCFU/mL
Chlamydia pneumoniae	≤1.0x10sCFU/mL
Legionella pneumophila	≤1.0x10sCFU/mL
Mycobacterium tuberculosis	≤1.0x10sCFU/mL
Pneumocystis yersinensis	≤1.0x10sCFU/mL
Pseudomonas aeruginosa	≤1.0x10sCFU/mL
Staphylococcus aureus	≤1.0x10sCFU/mL
Streptococcus salivarius	≤1.0x10sCFU/mL
Endemic human coronavirus HKU1	≤1 x 10s PFU/mL
SARS coronavirus	≤1 x 10s PFU/mL
Staphylococcus salivarius	≤1 x 10s CFU/mL
Human nasal rinse	≤14%

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## VIDEO TUTORIAL

## ANICE COVID-19 Antigen Saliva Lollipop (For Self-Test Only)

**REF** ISP-CA2 **Instruction for Use**  
**[Intended use]**

Used for the in vitro qualitative detection of the novel coronavirus (COVID-19) N protein antigen in human saliva samples. For Self-testing only. Antigen testing is generally used in the acute infection period, when samples are tested within seven days of the onset of symptoms in a suspected population. The positive result of the antigen test can be used for early triage and rapid management of suspected infected people, but the positive result only indicates the presence of the novel coronavirus N antigen in the sample, and cannot be used as the basis for the diagnosis and exclusion of pneumonia caused by the novel coronavirus. It should be combined with nucleic acid testing, imaging and other diagnostic information, medical history, and contact history to determine the status of infection. The negative result cannot rule out novel coronavirus infection, and it can not be used alone as a basis for treatment and disease management decisions. Suspected infected people with the antigen test results of positive and negative should be further nucleic acid testing. Coronavirus belongs to the order Nidoviridae, and the Coronavirus family is divided into three genera of  $\alpha$ ,  $\beta$ , and  $\gamma$ .  $\alpha$  and  $\beta$  are only pathogenic to mammals, and  $\gamma$  mainly causes bird infections.

CoV is mainly transmitted through direct contact with secretions or through aerosols and droplets. There is also evidence that it can be transmitted through fecal or oral routes. Until now, there are 7 kinds of human coronaviruses (HCoV) that cause human respiratory diseases: HCoV-229E, HCoV-OC43, SARS-CoV, HCoV-NL63, HCoV-HKU1, MERS-CoV and COVID-19, is an important pathogen of human respiratory infections. Among them, the clinical manifestations of SARS-CoV-2 virus are systemic symptoms such as fever and fatigue, accompanied by dry cough, dyspnea, etc., which can rapidly develop into severe pneumonia, respiratory failure, acute respiratory distress syndrome, septic shock, and multiple organ failure, severe acid-base metabolism disorder, etc., even life-threatening.

### 【Test principle】

Using double antibody sandwich method. The binding pad of this test strip was respectively covered with mouse anti-novel Coronavirus monoclonal antibody 1 with colloidal gold and Chicken IgY antibody as the color marker. The detection line (T) on the nitrocellulose membrane was covered with mouse anti-novel Coronavirus monoclonal antibody 2, and the quality control line (C) was covered with Goat anti-chicken IgY polyclonal antibody. When testing, when the sample to be tested contains the novel coronavirus, it combines with the colloidal gold-labeled novel coronavirus monoclonal antibody to form an immune complex, which is captured and enriched at the detection line (T) by the reagents fixed on the membrane. The colloidal gold-labeled antibody diffuses to the quality control line (C) area and is captured by the secondary antibody to form a purple-red band in the quality control area.

### 【Specifications and Main Components】

- 1 ANICE Covid-19 Antigen Saliva Lollipop
- 1 disposal biohazard bag and 1 disinfectant wipes
- 1 Instruction For Use

### Following materials required but not provided in the kit:

- Timer (watch)
- Any necessary personal protective equipment (gloves, glasses etc.)

### 【Storage and Shelf life】

The kit should be stored at 4–30°C. See the label for the production date and expiration date. After the reagent card is opened, it will not affect the performance of the kit within 60 minutes when exposed to temperature of 25°C and humidity of 40-65%.


### • 【Sample Collection】

1. Within 30 minutes before sampling, do not eating, smoking, drink wine and drink water. Wash the hands with hand sanitizer.
2. When taking a sample, the person should sit or stand firmly with the head level and slightly downward, gently put the absorbent tip of the test card into the mouth and place it under the tongue, so that the saliva is naturally absorbed by the absorbent tip.
3. During the sampling process, the test card (observation window must be upward), the level of the test card should be kept slightly downward, and the inclination angle of the test card should not be too high or too low (such as head down, head up).
4. Precautions: Saliva samples are suitable for testing this product, but it is recommended to collect saliva samples in the morning before rinsing, eating or drinking.
5. Samples heavily contaminated by oral food residues cannot be used for testing of this product. After a large amount of blood is contaminated, it is not recommended to use this product to test saliva samples. If the saliva sample is too viscous, the test results may vary greatly.

### 【Test Method】

Read the instructions carefully before using the kit, and operate in strict accordance with the instructions:

1. Take out the test reagent card, after equilibrating to room temperature (15-30°C), unpack the aluminum foil bag of the test reagent card and place it flat.
2. Remove the cover of the test card, put the adsorption head of the test card into the mouth, press the adsorption tip under the tongue, so that saliva is naturally absorbed by the tip of the adsorption tip, take out the test card after 2 minutes, and close the cover of the test card.
3. Place the test card horizontally, read the displayed result within 15-30 minutes, and the results read after 30 minutes is invalid

### 【Waste Disposal after Test Procedures】

1. After testing, place the used devices in the biohazard waste bag and seal the bag.

2. Throw away the remaining sample kit items.

3. Re-apply hand sanitizer.

### 【INTERPRETATION OF TEST RESULTS】


1. **Positive:** There are two red bands, and both the detection line (T line) and the quality control line (C line) are colored.

**Caution:** There is currently a suspicion of a COVID-19 infection. You are therefore encouraged to:

To contact a doctor/general practitioner or the local health department immediately. You need to self-report your result to My Sejahtera [https://mysejahtera.malaysia.gov.my/help\\_en/selfReport/](https://mysejahtera.malaysia.gov.my/help_en/selfReport/)

- Comply with local guidelines for self-isolation.
- Have a PCR confirmatory test performed.

2. **Negative:** Only a red band, the quality control line (C line) is colored;
3. **Invalid:** No color band appears on the quality control line (C), and it is judged as an invalid result regardless of whether the detection line (T) shows color band or not, and the test should be re-sampled.

**Positive Negative Invalid**

### 【Limitations】

1. This kit is a qualitative test and is only used for in vitro auxiliary diagnosis.
2. This kit is used for in vitro qualitative detection of novel coronavirus (COVID-19) N protein antigen in human saliva samples.
3. The positive test can only indicate that there may be a novel coronavirus (COVID-19) N protein antigen, and cannot be used as the only criterion for interpreting a novel coronavirus (COVID-19) infection. It should be combined with other diagnostic information such as nucleic acid testing and imaging Judging the infection status with medical history and contact history.
4. If the amount on N-protein antigen in a sample is below the limit of detection, a false-negative test result may occur.
5. Flawed sampling, transport, handling and insufficient virus content in a sample may lead to false-negative results.
6. Please ensure an appropriate sample size for testing. Too much or too little sample size may lead to inaccurate results.

### 【Performance Characteristics】

1. **Negative reference product coincidence rate**  
Take enterprise negative reference products (N1-N10) as samples for testing, and the coincidence rate should be 10/10.
2. **Positive reference product compliance rate**  
Take enterprise positive reference products (P1-P5) as samples for testing, and the coincidence rate should be 5/5.
3. **Reference Limit of Detection**  
S1~S4 are all positive, S5 can be positive or negative.
4. **Repeatability Reference**  
Take the corporate repeatability reference product (J) as the sample for parallel testing, and the test result should be positive.
5. **Cross reactivity:**  
The following pathogenic microorganisms were tested for cross reactivity with the COVID-19 Antigen Test Kit. For none of the tested microorganisms, cross reactivity was observed.

Analytes	Test concentration
Endemic human coronavirus NL63	≤9.87x10 <sup>6</sup> PFU/mL
MERS Coronavirus	≤793PFU/mL
Influenza A virus (H3N2)	≤8.82x10 <sup>4</sup> PFU/mL
Influenza B virus (Malaysia)	≤2.92x10 <sup>4</sup> PFU/mL
Rhinovirus	≤4.17x10 <sup>6</sup> PFU/mL
Endemic human coronavirus OC43	≤1.0x10 <sup>6</sup> PFU/mL
Endemic human coronavirus 229E	≤1.0x10 <sup>6</sup> PFU/mL
Adenovirus	≤1.0x10 <sup>6</sup> PFU/mL
Metapneumovirus	≤1.0x10 <sup>6</sup> PFU/mL
Parainfluenza virus type 1	≤1.0x10 <sup>6</sup> PFU/mL
Parainfluenza virus type 2	≤1.0x10 <sup>6</sup> PFU/mL
Parainfluenza virus type 3	≤1.0x10 <sup>6</sup> PFU/mL
Parainfluenza virus type 4	≤1.0x10 <sup>6</sup> PFU/mL
Influenza A virus (H1N1)	≤1.0x10 <sup>6</sup> PFU/mL
Enterovirus	≤1.0x10 <sup>6</sup> PFU/mL
Respiratory syncytial virus respectively	≤1.0x10 <sup>6</sup> PFU/mL
Haemophilus influenzae	≤1.0x10 <sup>6</sup> CFU/mL
Streptococcus pneumoniae	≤1.0x10 <sup>6</sup> CFU/mL
Streptococcus pyogenes	≤1.0x10 <sup>6</sup> CFU/mL
Candida albicans	≤1.0x10 <sup>6</sup> CFU/mL
Bordetella parapertussis	≤1.0x10 <sup>6</sup> CFU/mL
Mycoplasma pneumoniae	≤1.0x10 <sup>6</sup> CFU/mL
Chlamydia pneumoniae	≤1.0x10 <sup>6</sup> CFU/mL
Legionella pneumophila	≤1.0x10 <sup>6</sup> CFU/mL
Mycobacterium tuberculosis	≤1.0x10 <sup>6</sup> CFU/mL
Pneumocystis yersinensis	≤1.0x10 <sup>6</sup> CFU/mL
Pseudomonas aeruginosa	≤1.0x10 <sup>6</sup> CFU/mL
Staphylococcus aureus	≤1.0x10 <sup>6</sup> CFU/mL
Streptococcus salivarius	≤1.0x10 <sup>6</sup> CFU/mL
Endemic human coronavirus HKU1	≤1 x 10 <sup>6</sup> PFU/mL
SARS coronavirus	≤1 x 10 <sup>6</sup> PFU/mL
Staphylococcus salivarius	≤1 x 10 <sup>6</sup> CFU/mL
Human nasal rinse	≤14%

6. Interference : The following components were tested with COVID-19 Antigen Test Kit to determine their potential influence on the test outcome. All test results were negative, indicating that none of the components interfered with the test.

Analytes	Concentration
Purified mucin	≤10mg/mL
$\alpha$ -interferon	≤3million U
Ribavirin	≤2.0mg/mL
Oseltamivir	≤375ug/L
Peramivir	≤20ug/mL
Levofloxacin	≤5ug/mL
Azithromycin	≤0.15g/L
Ceftriaxone	≤100mg/mL
Meropenem	≤1ug/mL
Tobramycin	≤0.125mg/mL
Oxymetazoline	≤0.126ug/dL
Sodiumchloride	≤0.9%
Flutimethasone	≤200ug/mL
Budesonide	≤0.64nmol/L
Benzocaine	≤150mg/dL
Mupirocin	≤0.15mg/dL
Biotin	≤0.354mg/dL
Diphenhydramine	≤0.0774mg/dL
Dextromethorphan	≤0.00156mg/dL
Dexamethasone	≤1.2mg/dL
Normal nasopharyngeal mucus	≤5%
Whole blood	≤5%

7. Limit of Detection: Use normal human saliva samples as the negative matrix diluent to detect the novel coronavirus inactivated virus solution after serial dilution. The lowest detection limit is 600 TCID50/mL.
8. Hook effect: The samples in the test that the concentration level of the novel coronavirus inactivated virus solution is 2.0x10<sup>5</sup>TCID50/mL did not show the Hook effect.

### 【Clinical Performance】

The results of this product and the novel coronavirus nucleic acid detection reagent (RT-PCR) were compared and studied. 821 saliva samples were tested. The results are shown below.

Test	RT-PCR		Total
	Positive	Negative	
ANICE COVID-19 Antigen Saliva Lollipop (For Self-Test Only)	390	3	393
	15	413	428
Total	405	416	821

Positive coincidence rate (Sensitivity): 96.30% (95% CI: 93.96%~97.91%)  
 Negative coincidence rate (Specificity): 99.28%(95%CI: 97.91%~99.85%)  
 Total coincidence rate: 97.81% (95%CI: 96.56%~98.70%)

### 【Precautions】

1. This kit is only used for in vitro testing. Please read this manual carefully before experimenting and strictly follow the operating procedures in the manual.
2. The collection, storage and testing of samples should be carried out in strict accordance with relevant guidelines.
3. After the inspection, the remaining samples are preserved and various waste treatments. The waste or remaining samples generated during the inspection process are recommended to refer to the above guidelines. First, diethyl ether, 75% ethanol, chlorine-containing disinfectant, peracetic acid, chloroform and other lipid solvents are used to soak the virus for inactivating, and then the infectious agents are treated according to the above guidelines.

### 【Frequently Asked Questions (FAQ)】

- **When can/should I test myself?**
- You can have a test on yourself whether you have symptoms or not. Please note that the test result is a snapshot that is valid for this point in time. Tests should therefore be repeated according to local regulations.
- **I have taken the test, but the control line (C) doesn't appear. What should I do?**
- According to the instructions for use, this test result is invalid. Please retest by using a new test card.
- **I am not sure about reading test result. What should I do?**
- Read the instructions for use again, and if this doesn't help, please contact the nearest health facility recommended by your local authorities for help.
- **If my test result is positive, what should I do?**
- There is possibility of hospitalization, complications and even death after infection with COVID-19. You should immediately contact the nearest health facility recommended by your local authorities.
- **If my test result is negative, what should I do?**
- If you test result is negative by the test, you also need to obey the local regulations. If you experience symptoms such as fever, headaches, migraines, loss of sense of smell and taste, contact the nearest health facility recommended by your local authorities.

### 【Symbol】

Symbol	Explanation	Symbol	Explanation
	In vitro diagnostic medical device		Lot number
	Contains sufficient for <n>tests		Date of manufacture
	Manufacturer		Expiry Date
	Authorized representative in the European Community		Temperature limit: 4~30°C
	CE Marking		This way up
	Caution		Keep dry
	Consult instructions for use		Fragile: handle with care
	Biological risks		Keep away from sunlight
	Catalogue number		Do not re-use

ISP WELLNESS SDN BHD (1433750-M)

No 47, Jalan Chung Ah Ming, 31650 Ipoh, Perak, Malaysia.

Tel: +(60)5-321 7673

Thank you for purchasing ANICE COVID-19 Antigen Saliva Lollipop (for self-test only). Please read this manual carefully before operating to ensure proper use.

VIDEO TUTORIAL